

Rehabilitation Treadmill - 7.0T



Multiple adjustments
for individualized fit

The Rehabilitation Treadmill - 7.0T is packed with features that reflect the true needs of physical therapists, this translates into patients feeling confident and comfortable in using this for their recovery in physical mobility.

Benefits



The step makes access to the deck easier

The removable step provides a lower step-up height for people who need a bit more stability getting on and off the treadmill. This is ideal for individuals who cannot make a big step - such as individuals with balance limitations, smaller range of motion, or older adults.



Multiple default and customisable programs are available

With many features, this easy-to-use console displays and downloads important workout data. Feedback windows display METs, symmetry, cadence, stride lengths, steps, and pace. You can easily select either preset or programmable modes for diverse, versatile routines. Programs can interactively adapt to patients' characteristics, using input info such as weight and gender as the basis to maximise their rehabilitation.



Fully adjustable parallel bars supports a wide variety of users

The adjustable parallel bars offer handrail support for added patient safety. Adjustments can also accommodate users of various heights and fitness abilities. Bars adjust 63.5 to 89 cms in height, 51cm to 85cm in width, and the length extends to the full length of the deck.

Specifications

Incline	-10% to 15%
Bidirectional Speed	Forward 0.1 to 16 kph / 10 mph and reverse 0.1 to 4.8 kph / 3 mph in 0.1 increments
Motor	3HP continuous duty, IEC Class S6
Belt and Deck	152.5 x 56 cm / 22" x 60" self-lubricating deck
Power	115 volts AC, 60Hz, 20 amps, NEMA-20P
Max User Weight	200 kg / 440 lbs
Dimensions	213 x 144 x 140 cm / 83.7" x 56.7" x 55"
Product Weight	174 kg / 384 lbs
Compliance	CAN/CAS-C22.2 No. 60601-1:14 ANS/AAMI ES60601-1:2005+A2 (R2012) +A1 IEC 60601-1-2:2014 EN 60601-1-2:2015 IEC 60601-1:2015+A1:2012 EN 60601-1:2006+A1:2013+A12:2014 EN ISO 20957:2013 MDD 93/42/EEC Class IIa, NB: 0123
Warranty	Provided upon request with your local distributor

For more information contact your local distributor

<https://medical.dyaco.com/distributors>